DASHBAORD PRESENTATION





Presented by: D. pavan, R.kavya

Key Functionalities

- Store JSON Data in MongoDB
- Upload the pdf files from the frontend
- Generate Excel Files
- Submit the reports
- Send Emails with Attachments

Flowchart

_id: ObjectId('67c97fc16d4743d1fd335336')

batch_name: "BATCH1"

▼ patients: Array (6)

▼ 0: Object

patient_id : "KHAIGHGPPGX27"

▶ files: Object

▼ 1: Object

patient_id : "KHAIGHGPTTL569"

▶ files: Object

▶ 2: Object
▶ 3: Object

▶ 4: Object

▶ 5: Object



Pathogenic Variants

Pathogenic Variants

rsID ↑↓	Gene Name ↑↓	Gene Score ↑↓	Lit ↑↓	сн↑↓	pos ↑↓	ref ↑↓	alt ↑↓	Zygosity ↑↓	Consequence ↑↓	Conseq ↑↓	clin ↑↓	IMPACT ↑↓
rs543016186	MTHFR	No	No	chr1	11795125	С	Т	Heterozygous	missense variant	7/10	Pathogenic	MODERATE
rs201420507	DIO1	No	No	chr1	53906216	G	Α	Heterozygous	missense variant	7/10	Pathogenic	MODERATE
rs1559416138	WNT10A	No	No	chr2	218892949	Α	-	Heterozygous	frameshift variant	10/10	Pathogenic	HIGH
rs3974499	SLC9B1	No	No	chr4	102901326	AC	1-1	Heterozygous	frameshift variant	10/10	Pathogenic	HIGH
rs199422117	TBXAS1	No	No	chr7	140015731	G	Α	Heterozygous	missense variant	7/10	Pathogenic	MODERATE
rs202003805	PRSS1	8	Yes	chr7	142750561	С	T	Heterozygous	missense variant	7/10	Pathogenic	MODERATE
rs786204713	ACADVL	7	No	chr17	7220783	AC	-	Heterozygous	frameshift variant	10/10	Pathogenic	HIGH
rs60726084	ATXN10	No	No	chr22	45795355	ATTCT	(=)	Heterozygous	intron variant	2/10	Pathogenic	MODIFIER


Dashboard Overview Content:

• **Purpose:** The dashboard provides a user-friendly interface to manage patient reports and their availability status.

• Key Features:

- Submit patient reports
- Al Score Display
- Severity Selection
- Full screen toggle mode
- Update patient availability
- View report submission and availability status
- Upload the pdf files directly into the database from frontend

Dashboard



DashboardSearch Component

1. Batch Selection:

- Dropdown to select a batch.
- Fetches batches from API on mount.
- Updates UI and fetches report status when a batch is selected.

2. Patient Search:

- Filters patients dynamically as the user types.
- Displays all patients if search input is empty or focused.
- Allows selecting a patient to view details.

3. Patient Availability Toggle:

- Toggles patient availability via a switch.
- Prevents toggling if the report is already submitted.
- Sends updates to the backend and updates local state.


4. Report Status Handling:


- Fetches and displays report submission status.
- Uses color coding:
 - Green: Report submitted.
 - Orange: Available.
 - Red: Not available.

5. View Details Button:

- Allows users to select a patient and view details.
- Calls onSelectPatient and fetchDataFromAPIs upon selection.


Patient selection, availability status & Report submisssion indication





Condition button colour indication

- If the condition button was coloured in the orange colour it indicates that condition was having concern.
- If the condtion button was coloured in light green colour it indicates currently that condition was selected.
- If the condition button was coloured in the Dark green colour it indicates that condition was viewed and given report.



Full-screen mode


Key Features:

- Full-Screen Toggle: Enter/exit fullscreen mode.
- Severity Selection: Buttons with icons & color-coded levels.
- Concern Checkbox: Mark/unmark a condition as a concern.
- AI Score Display: Toggle AI-generated score visibility.
- No Mutation Button: Indicates absence of mutation.
- Condition Details: Dynamically renders selected conditions.
- Responsive Design: Scrollable and adaptable UI.

Full-screen mode


		②	Low		ld 😉) Moderate	⊗ v	Moderate to High 🔽 🤇	Concern Al Score	X No Mutatio	on	
Pancreatitis												
Pancreatitis												
rsID ↑↓	Gene Name ↑↓	Gene Score ↑↓	Lit ↑↓	сн ↑↓	pos ↑↓	ref ↑↓	alt 👈	Zygosity ↑↓	Consequence ↑↓	Conseq ↑↓	clin ↑↓	IMPACT 1
rs202003805	PRSS1	8	Yes	chr7	142750561	С	Т	Heterozygous	missense variant	7/10	Pathogenic	MODERATE
rs144422014	PRSS1	8	Yes	chr7	142750675	Α	G	Heterozygous	missense variant	7/10	Likely Benign	MODERATE
rs213950	CFTR	8	No	chr7	117559479	G	Α	Homozygous Variant	missense variant	7/10	Likely Benign	MODERATE
rs748442280	PRSS1	8	No	chr7	142749516	С	G	Heterozygous	missense variant	7/10	Likely Benign	MODERATE
rs747228052	PRSS1	8	No	chr7	142749524	С	G	Heterozygous	missense variant	7/10	Conflicting	MODERATE
rs200665515	PRSS1	8	No	chr7	142750558	С	G	Heterozygous	missense variant	7/10		MODERATI
rs770782578	PRSS1	8	No	chr7	142750563	С	Т	Heterozygous	missense variant	7/10	Likely Benign	MODERATI
rs138464021	PRSS1	8	No	chr7	142750660	G	Т	Heterozygous	missense variant	7/10	Benign	MODERATE
rs149246646	PRSS1	8	No	chr7	142750672	T	Α	Heterozygous	missense variant	7/10		MODERATI
rs148440491	PRSS1	8	No	chr7	142750676	С	G	Heterozygous	missense variant	7/10		MODERATI
rs147366981	PRSS1	8	No	chr7	142750680	С	T	Heterozygous	stop gained	10/10	Conflicting	HIGH
rs757111793	PRSS1	8	No	chr7	142751776	G	A	Heterozygous	missense variant	7/10	Uncertain Significance	MODERATE
rs750348889	PRSS1	8	No	chr7	142751865	С	Α	Heterozygous	missense variant	7/10	Likely Benign	MODERATI
rs199507985	PRSS1	8	No	chr7	142751871	G	Α	Heterozygous	missense variant	7/10	Conflicting	MODERATI
rs778568523	PRSS1	8	No	chr7	142751878	Α	Т	Heterozygous	missense variant	7/10	Uncertain Significance	MODERATI
rs1373631104	PRSS1	8	No	chr7	142751879	G	Т	Heterozygous	missense variant	7/10	Uncertain Significance	MODERATI
rs1232891794	PRSS1	8	No	chr7	142752476	G	С	Heterozygous	missense variant	7/10	Uncertain Significance	MODERATI
rs758254763	PRSS1	8	No	chr7	142752490	G	Α	Heterozygous	missense variant	7/10		MODERATI
rs1240508430	PRSS1	8	No	chr7	142752505	G	Т	Heterozygous	stop gained	10/10	Uncertain Significance	HIGH
rs1468060476	PRSS1	8	No	chr7	142752506	G	Т	Heterozygous	missense variant	7/10	Uncertain Significance	MODERATE
rs1348773645	PRSS1	8	No	chr7	142752522	С	G	Heterozygous	missense variant	7/10	Uncertain Significance	MODERATE
rs200902389	PRSS1	8	No	chr7	142752913	G	Α	Heterozygous	missense variant	7/10	Benign	MODERATE

PDF Upload



 Added one file upload button in the navbar that allows the user to upload the pdf reports like Blood-work, Consent, directly into the database.

Multiconditions Selection



 Added one Multicondition button in the navabar which allows the user to view the two condtions one patient at a time.

PDF Viewer

PDF Display:

- It loads and displays a multi-page PDF document.
- The current page is "1 / 44," indicating a total of 44 pages.


Navigation Controls:


- Buttons to move between pages (+ and zoom controls).
- A refresh/reset button to reload the document.

Download & Options:

- A download button (down arrow icon) to save the PDF.
- A settings/menu button (three dots icon) for more options.
- Checkboxes & Buttons:
- Buttons labeled PDF, Consent, Blood Reports to categorize documents.
- Scrollable View:
- A scroll bar on the right allows vertical navigation.

PDF Viewer






<u>Submitting Patient Reports Content:</u>

• Functionality:

- Users can submit patient reports by filling out a form with necessary details such as patient ID, batch, and report data.
- Upon submission, the data is stored in MongoDB.
- An Excel file is generated from the submitted data.
- The generated Excel file is sent via email to specified recipients.
- The file is deleted from the server after the email is sent.

Submitting Patient Report



Submitted Reports in Email

X KHAIGHGPPGX27_Scoring_chart.xlsx												Open w
	А	В	С	D	E	F	G	Н	ı	J	К	L
1	dical Conditi	Low	Low to Mild	Mild	ld to Modera	Moderate	oderate to Hi	High	concerns	No Mutations	Al Score	Reason
2	Cardiac_Hea	у							у	у	Moderate to	nothingj
3	Cholesterol_	_Disorders				у			у		Mild	nothingj
4	Diabetes	у							у	у	Moderate	test
5	High_Blood_	Pressure				у			у		Low	last
6	Gastritis			у					у		Moderate	final
7	Glomerular_	у							у	у	Mild	nothing
8	Gut_Health	у							у		Moderate to	final
9	Headaches	у							у		Moderate	normal t
10												
11												
12												
13												
14												

LIVE DEMO

